

FDA Submission


Gaurav Yadav

Pneumonia screener

Algorithm Description

1. General Information

The algorithm was trained on 56.5% male and 43.5 female patients whom fall in range of 5 to 90. All the patient were scanned for the **chest X-ray** and were labeled with 14 disease and **No Finding**. From the below plot, it is clear that most of the patients are from age group 50–65. And we can say distribution among gender is almost balance.



Intended Use Statement:

The algorithm is intended to use for Pneumonia patient who has been adminstered to a chest X-ray screening and had never demonstrated a chest abnormalities.

Indications for Use:

The algorithm can be used for screening the pneumonia patient using chest x-ray which can be helpful in early detection of pneumonia. This algorithm would be helpful to the radiologist to filter out the patient who doesn't have pneumonia.


Device Limitations:

The algorithm is to detect whether patient has pneumonia or not but this can perform poorly in presence of other disease like Atelectasis, Cardiomegaly, Consolidation, Edema, Effusion, Emphysema, Fibrosis, Hernia, Infiltration, Mass, Nodule, Pleural Thickening, or Pneumothorax. If the algorithm say patient has pneumonia then it is possible that it may not have pneumonia as we are focussing more on the recall instead of precision.

The algorithm is trained using chest x-ray images with view angle AP or PA. So, it's optimal performance would be while using chest x-ray image which are take with view angle of either PA or AP. It will perform poorly on other inputs like CT etc.

Clinical Impact of Performance: The algorithm can be great help to the radiologist for studying chest xray of a patient. As the algorithm could filter out the patient which don't have pneumonia. Radiologist will have to look at less number of chest xrays.

2. Algorithm Design and Function

Basic flow chart of the system


DICOM Checking Steps: Below is a list of DICOM checks:

1. DICOM file must be of **chest** body part. (`BodyPartExamined == 'CHEST' .`)
2. DICOM file should cotain XRAY image only. (`Modality == 'DX' .`)
3. DICOM image view angle must be either PA or AP. (`PatientPosition in ('PA', 'AP') .`)

Preprocessing Steps: Training data contains images with size (1024,1024). For training data has been preprocessed as follow:

- Image has been resized to (224,224)
- Image has been auggedmented wiht the following parameters
 - Image rescaled to 1/255 to have pixel value in range of 0-1
 - Rotation range of 20
 - Zoom range of 0.1
 - Shear range of 0.1
 - Horizontal flip applied

CNN Architecture: The algorithm has been designed using VGG16 pretrained model which was initially developed for imagenet challenge. Below is the architectural image of VGG and after that Algorithm's structured has been appended. After last convolution layer of VGG16, one flatten layer has been added then a batchNormalization then a Dense layer folloing a Dropout layer and lastly Dense layer with single neuron.




Final model:

Layer (type)	Output Shape	Param #
model_2 (Model)	(None, 7, 7, 512)	14714688
flatten_2 (Flatten)	(None, 25088)	0
batch_normalization_2 (Batch Normalization)	(None, 25088)	100352
dense_3 (Dense)	(None, 256)	6422784
dropout_2 (Dropout)	(None, 256)	0
dense_4 (Dense)	(None, 1)	257


3. Algorithm Training

Parameters:

- Types of augmentation used during training
- Batch size = 100
- Optimizer learning rate = $1e^{-4}$
- Layers of pre-existing architecture that were frozen
 - first 17 layers of VGG has been frozen
- Layers of pre-existing architecture that were fine-tuned
 - block5_pool
- Layers added to pre-existing architecture
 - Flatten
 - BatchNormalization
 - Dense with activation relu
 - Dropout with probability 0.5
 - Dense (last layer with 1 neuron) with sigmoid activation.



And the precision-recall curve has an AP score of 0.47.



The algorithm has an area under the curve for True positive rate and false positive rate of 0.70.

Final Threshold and Explanation:

Here, we are specially focused on reduction the number of False positive (FP) and False Negatives(FN), mostly False Negative, of confusion matrix as the numbers of cases under False Positives would conclude that the patients, who are actually Normal, are facing Pneumonia, whereas False Negatives would conclude that the patients, who actually have Pneumonia, are not classified as suffering from Pneumonia. **Since we can tolerate the number of cases of False Positives, but not False Negatives! This may have serious repercussions if not classified correctly. So it is of utmost importance to reduce the number of cases of False Negatives.**

So, our focus is to increase the recall instead of precision. Below is the threshold and f1 score for recall of 0.8

```
Precision is: 0.31666666666666665
Recall is: 0.7916666666666666
Threshold is: 0.22881234
F1 Score is: 0.45238095238095233
```

4. Databases

This NIH Chest X-ray Dataset is comprised of 112,120 X-ray images with disease labels from 30,805 unique patients. To create these labels, the authors used Natural Language Processing to text-mine disease classifications from the associated radiological reports. The labels are expected to be >90% accurate and suitable for weakly-supervised learning. The original radiology reports are not publicly available but you can find more details on the labeling process in this Open Access paper: "ChestX-ray8: Hospital-scale Chest X-ray Database and Benchmarks on Weakly-Supervised Classification and Localization of Common Thorax Diseases." (Wang et al.)

[Link to paper](#)

Data limitations:

- The image labels are NLP extracted so there could be some erroneous labels but the NLP labeling accuracy is estimated to be >90%.
- Very limited numbers of disease region bounding boxes.
- Data is highly imbalanced and has very few pneumonia positive data point. We only have 1.2% of pneumonia patient in whole dataset i.e. 1430/112104

Labels and patient data for the entire dataset

- Image Index: File name
- Finding Labels: Disease type (Class label)
- Follow-up #
- Patient ID
- Patient Age
- Patient Gender
- View Position: X-ray orientation
- OriginalImageWidth
- OriginalImageHeight
- OriginalImagePixelSpacing_x
- OriginalImagePixelSpacing_y

Class descriptions


There are 15 classes (14 diseases, and one for "No findings"). Images can be classified as "No findings" or one or more disease classes:

- Atelectasis
- Consolidation
- Infiltration
- Pneumothorax
- Edema
- Emphysema
- Fibrosis
- Effusion
- Pneumonia
- Pleural_thickening
- Cardiomegaly
- Nodule Mass
- Hernia

Description of Training Dataset: Below are some properties from the training data.


Gender Disgribution

Gender distribution in dataset




Age Distribution


Hist. Age distribution	Age bins distribution
------------------------	-----------------------




View angle Distribution



Disease Distribution



Pneumonia Co-occurrence with other disease



It is quite clear from the above analysis, data set is highly imbalanced and has only few positive pneumonia patients. But for the training purpose we are not using the whole data set but a balanced data where number of pneumonia patients is equal to the non-pneumonia patients. Also some outlier data points like age > 110 have been dropped. And all the images have been preprocessed before being used by the algorithm.

Description of Validation Dataset: Dataset has been divided into train and validation sets following the 80-20 ratio. The validation data set doesn't have the same image preprocessing of training. In validation set, we are only resizing the image and normalizing the pixel value, no other additional pre-processing has been applied.

5. Ground Truth

The biggest limitation of this dataset is that image labels were NLP-extracted so there could be some erroneous labels but the NLP labeling accuracy is estimated to be >90%. We can have optimal ground truth if radiologists examine the samples and label them accordingly.

6. FDA Validation Plan

Patient Population Description for FDA Validation Dataset: For the validation of this algorithm, data should be collected from men and women distributed between the age range of 1-110 with no prior history of pneumonia and other chest abnormal diseases like Edema, Effusion, Emphysema, Fibrosis, Hernia, Infiltration, Mass, Nodule, Pleural, Atelectasis, Cardiomegaly, Consolidation, Thickening, or Pneumothorax.

Data should be collected as DICOM files with the following properties:

1. BodyPartExamined must be chest as we have trained the algorithm on chest images only.
2. Image view angle must be either AP or PA. In DICOM file, PatientPosition must be either AP or PA.
3. Images are of XRAY type only, means Modality of DICOM file must be 'DX'

Ground Truth Acquisition Methodology:

For ground truth, we need to label each of our data collected with class label whether the x-ray image is of pneumonia patient or Non-Pneumonia (Healthy) patient. This can be achieved following the silver standard approach. So, several radiologists can label the collected data with class label (Pneumonia or Healthy (non-pneumonia)). For more optimal ground truth, we should also consider the experience of radiologists.

Algorithm Performance Standard:

Below is the comparison of our algorithm with others in terms of AUROC:

Pathology	Wang et al. (2017)	Yao et al. (2017)	CheXNet	ours
Pneumonia	0.633	0.713	0.7680	0.71

Ref: [paper](#)